

DNA Isolation from Saliva FTA

DNA extraction from saliva FTA using *forensicGEM*[™]

The following method is recommended for DNA extraction from FTA cards using *forensicGEM*[™].

The method below is intended for a 96 well format. However, the same method is applicable for any number of samples using a PCR tube and a thermal cycler.

Extraction Method


1. Add 1 x 1.2 mm FTA punch to each well of a microtitre dish [NOTE1].
2. Add 50 μ l of *forensicGEM*[™] buffer 3 and 1 μ l of *forensicGEM*[™] [NOTE1].
3. Incubate at 75 °C for 15 minutes.
4. Incubate at 95°C for 5 minutes [NOTE2].
5. Quantify the supernatant using Quantiblot or Quantifiler if required.

NOTE1: Larger punch sizes can be used with the same method. The yield will be proportionately higher. *forensicGEM*[™] and buffer can also be added as a master mix.

NOTE2: 5 minutes is sufficient as long as ramping to temperature is rapid. For slow ramping, longer incubations should be used.

SALIVA FTA

The results below were from a PCR using 5 μ l of an extraction derived from a single 1.2 mm FTA card punch of saliva swabs. The top panel compares profiles using AmpFLSTR[®] Identifiler[™] PCR Amplification Kit (Applied Biosystems). The bottom panel shows results using SGM+.



Comparison of yields

DNA extractions with *forensicGEM*[™] produced slightly lower yields overall but all gave comparable profiles.

