

Easy-Go™ RT PreMix

PROTOCOL:

1. Mix the template RNA and primer in a sterile tube as indicated below:

Reaction volume		20 µl reaction
<i>Template RNA</i>	Total RNA	0.5-1.0 µg
	Poly(A) ⁺ RNA	0.05-0.1 µg
<i>Primer</i>	Sequence specific	10-30 pmol
	Oligo dT ₁₈	100 pmol
	Random primers	100 pmol

2. Incubate the mixture at **70°C for 5 min** and place it on ice.
3. Transfer the incubated mixture to an **Easy-Go™ RT PreMix tube**, then make up the reaction volume with DEPC-Water to 20 µl.
4. Dissolve the lyophilized pellet by vortexing, and briefly spin down.
5. Add mineral oil to each tube (this step is unnecessary when using a thermal cycler with top heating).
6. Perform cDNA synthesis reaction as follows:
42 °C, 60 min (cDNA synthesis)
94 °C, 5 min.(RTase inactivation)

If PCR is needed following RT reaction, perform the PCR reaction as follows:

1. Transfer 2-5 µl of the RT product (synthesized cDNA) to **Easy-Do™ PCR PreMix tube** (20 µl reaction volume).
2. Perform PCR cycles according to the PCR condition.